

Cma1 Antibody, FITC conjugated



PACO25822

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Major secreted protease of mast cells with suspected roles in vasoactive peptide generation, extracellular matrix degradation, and regulation of gland secretion.

Gene ID:

Cma1

Uniprot

P21844

Synonyms:

Chymase (EC 3.4.21.39) (Alpha-chymase) (Mast cell chymase 1) (Mast cell protease 5) (mMCP-5) (Mast cell protease I), Cma1, Mcpt5

Immunogen:

Recombinant Mouse Chymase protein (22-246AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A