

RPL14 Antibody, HRP conjugated



PACO25838

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins.

Gene ID:

RPL14

Uniprot

P50914

Synonyms:

60S ribosomal protein L14 (CAG-ISL 7) (Large ribosomal subunit protein eL14), RPL14

Immunogen:

Recombinant Human 60S ribosomal protein L14 protein (2-214AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A