

# RPL14 Antibody, FITC conjugated



PACO25839

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins.

**Gene ID:**

RPL14

**Uniprot**

P50914

**Synonyms:**

60S ribosomal protein L14 (CAG-ISL 7) (Large ribosomal subunit protein eL14), RPL14

**Immunogen:**

Recombinant Human 60S ribosomal protein L14 protein (2-214AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A