

COL11A1 Antibody, HRP conjugated



PACO25850

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play an important role in fibrillogenesis by controlling lateral growth of collagen II fibrils.

Gene ID:

COL11A1

Uniprot

P12107

Synonyms:

Collagen alpha-1(XI) chain, COL11A1, COLL6

Immunogen:

Recombinant Human Collagen alpha-1(XI) chain protein (532-699AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A