

# CPA2 Antibody, FITC conjugated



PACO25895

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

extracellular region, extracellular space, vacuole, carboxypeptidase activity, metalloproteinase activity, zinc ion binding, protein catabolic process in the vacuole, proteolysis.

**Gene ID:**

CPA2

**Uniprot**

P48052

**Synonyms:**

Carboxypeptidase A2 (EC 3.4.17.15), CPA2

**Immunogen:**

Recombinant Human Carboxypeptidase A2 protein (115-419AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A