

PACO25905

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Receptor for complement C3d. Participates in B lymphocytes activation.

**Gene ID:**

Cr2

**Uniprot**

P19070

**Synonyms:**

Complement receptor type 2 (Cr2) (Complement C3d receptor) (CD antigen CD21), Cr2

**Immunogen:**

Recombinant Mouse Complement receptor type 2 protein (729-963AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A