

CRABP2 Antibody, FITC conjugated



PACO25915

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transports retinoic acid, to the nucleus. Regulates the access of retinoic acid, to the nuclear retinoic acid, receptors.

Gene ID:

CRABP2

Uniprot

P29373

Synonyms:

Cellular retinoic acid, binding protein 2 (Cellular retinoic acid, binding protein II) (CRABP-II), CRABP2

Immunogen:

Recombinant Human Cellular retinoic acid, binding protein 2 protein (1-138AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A