

# CRADD Antibody, HRP conjugated



PACO25918

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Apoptotic adaptor molecule specific for caspase-2 and FASL/TNF receptor-interacting protein RIP. In the presence of RIP and TRADD, CRADD recruits caspase-2 to the TNFR-1 signalling complex.

**Gene ID:**

CRADD

**Uniprot**

P78560

**Synonyms:**

Death domain-containing protein CRADD (Caspase and RIP adapter with death domain) (RIP-associated protein with a death domain), CRADD, RAIDD

**Immunogen:**

Recombinant Human Death domain-containing protein CRADD protein (1-199AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A