

CS Antibody, Biotin conjugated



PACO25968

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Citrate synthase is a pace-maker enzyme in the Krebs cycle.

Gene ID:

CS

Uniprot

O75390

Synonyms:

Citrate synthase, mitochondrial (EC 2.3.3.1) (Citrate (Si)-synthase), CS

Immunogen:

Recombinant Human Citrate synthase, mitochondrial protein (28-464AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A