

# CSNK2A1 Antibody, FITC conjugated



PACO26003

---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	CSNK2A1
<b>Source:</b>	<b>Uniprot</b>
Rabbit	P68400
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	Casein kinase II subunit alpha (CK II alpha) (EC 2.7.11.1), CSNK2A1, CK2A1
<b>Applications:</b>	<b>Immunogen:</b>
ELISA	Recombinant Human Casein kinase II subunit &alpha; protein (1-391AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A