

PACO26029

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

This vimentin-specific protease may regulate the reorganization of vimentin filaments, occurring during cell differentiation, movement and mitosis.

Gene ID:

Ctsg

Uniprot

P28293

Synonyms:

Cathepsin G (EC 3.4.21.20) (Vimentin-specific protease) (VSP), Ctsg

Immunogen:

Recombinant Mouse Cathepsin G protein (21-261AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A