

CTSL Antibody, Biotin conjugated



PACO26040

Product Information

Size:

50ug

Reactivity:

Chlorocebus aethiops

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Important for the overall degradation of proteins in lysosomes.

Gene ID:

CTSL

Uniprot

Q9GKL8

Synonyms:

Cathepsin L1 (EC 3.4.22.15) (Cathepsin L) (Major excreted protein) (MEP) [Cleaved into: Cathepsin L1 heavy chain; Cathepsin L1 light chain], CTSL, CTSL1

Immunogen:

Recombinant Chlorocebus aethiops Cathepsin L1 protein (114-333AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A