

BRINP1 Antibody, FITC conjugated



PACO26143

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Inhibits cell proliferation by negative regulation of the G1/S transition. Mediates cell death which is not of the classical apoptotic type and regulates expression of components of the plasminogen pathway.

Gene ID:

BRINP1

Uniprot

O60477

Synonyms:

BMP/retinoic acid, inducible neural-specific protein 1 (Deleted in bladder cancer protein 1), BRINP1, DBC1, DBCCR1, FAM5A

Immunogen:

Recombinant Human BMP/retinoic acid, inducible neural-specific protein 1 protein (20-320AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A