

DCN Antibody



PACO26153

Product Information

Size:

50ug

Reactivity:

Bovine, Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

May affect the rate of fibrils formation.

Gene ID:

DCN

Uniprot

P21793

Synonyms:

Decorin (Bone proteoglycan II) (PG-S2), DCN


Immunogen:

Recombinant Bovine Decorin protein (31-360AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: DCN antibody at 6 μ g/ml + Jurkat whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 40 kDa

Observed band size: 40 kDa