

RPS19 Antibody



PACO26161

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

Required for pre-rRNA processing and maturation of 40S ribosomal subunits.

Gene ID:

RPS19

Uniprot

P39019

Synonyms:

40S ribosomal protein S19 (Small ribosomal subunit protein eS19), RPS19


Immunogen:

Recombinant Human 40S ribosomal protein S19 protein (1-145AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot
All lanes: RPS19 antibody at 2 μ g/ml
Lane 1: EC109 whole cell lysate
Lane 2: 293T whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/15000 dilution
Predicted band size: 17 kDa
Observed band size: 17 kDa
.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO26161 at dilution of 1: 50.