

# DDAH2 Antibody, Biotin conjugated



PACO26180

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Hydrolyzes N(G), N(G)-dimethyl-L-arginine (ADMA) and N(G)-monomethyl-L-arginine (MMA) which act as inhibitors of NOS. Has therefore a role in the regulation of nitric oxide generation.

**Gene ID:**

DDAH2

**Uniprot**

O95865

**Synonyms:**

N(G), N(G)-dimethylarginine dimethylaminohydrolase 2 (DDAH-2) (Dimethylarginine dimethylaminohydrolase 2) (EC 3.5.3.18) (DDAHII) (Dimethylargininase-2) (Protein G6a) (S-phase protein), DDAH2, DDAH G6A NG30

**Immunogen:**

Recombinant Human N(G), N(G)-dimethylarginine dimethylaminohydrolase 2 protein (1-285AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A