

Ddt Antibody, HRP conjugated



PACO26190

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Tautomerization of D-dopachrome with decarboxylation to give 5,6-dihydroxyindole (DHI).

Gene ID:

Ddt

Uniprot

O35215

Synonyms:

D-dopachrome decarboxylase (EC 4.1.1.84) (D-dopachrome tautomerase), Ddt

Immunogen:

Recombinant Mouse D-dopachrome decarboxylase protein (2-118AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A