

# Ddt Antibody, FITC conjugated



PACO26191

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Tautomerization of D-dopachrome with decarboxylation to give 5,6-dihydroxyindole (DHI).

**Gene ID:**

Ddt

**Uniprot**

O35215

**Synonyms:**

D-dopachrome decarboxylase (EC 4.1.1.84) (D-dopachrome tautomerase), Ddt

**Immunogen:**

Recombinant Mouse D-dopachrome decarboxylase protein (2-118AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A