

DDX58 Antibody, Biotin conjugated



PACO26200

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

DDX58

Uniprot

O95786

Synonyms:

Probable ATP-dependent RNA helicase DDX58 (EC 3.6.4.13) (DEAD box protein 58) (RIG-I-like receptor 1) (RLR-1) (Retinoic acid, inducible gene 1 protein) (RIG-1) (Retinoic acid, inducible gene 1 protein) (RIG-I), DDX58

Immunogen:

Recombinant Human Probable ATP-dependent RNA helicase DDX58 protein (1-430AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A