

Defa5 Antibody, Biotin conjugated



PACO26204

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probably contributes to the antimicrobial barrier function of the small bowel mucosa.

Gene ID:

Defa5

Uniprot

P28312

Synonyms:

Alpha-defensin 5 (Defensin-related cryptdin-5), Defa5, Defcr5

Immunogen:

Recombinant Mouse Alpha-defensin 5 protein (59-93AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A