

# DEFB126 Antibody



PACO26209

---

## Product Information

**Size:**

50ug

**Reactivity:**

Pongo pygmaeus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Has antibacterial activity.

**Gene ID:**

DEFB126

**Uniprot**

A4H244

**Synonyms:**

Beta-defensin 126 (Defensin, beta 126), DEFB126

**Immunogen:**

Recombinant Pongo pygmaeus Beta-defensin 126 protein (21-63AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A