

GSDME Antibody, FITC conjugated



PACO26223

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in apoptosis and cell survival. Plays a role in the TP53-regulated cellular response to DNA damage probably by cooperating with TP53.

Gene ID:

GSDME

Uniprot

O60443

Synonyms:

Gasdermin-E (Inversely correlated with estrogen receptor expression 1) (ICERE-1) (Non-syndromic hearing impairment protein 5) [Cleaved into: Gasdermin-E, N-terminal (GSDME-NT); Gasdermin-E, C-terminal (GSDME-CT)], GSDME, DFNA5 ICERE1

Immunogen:

Recombinant Human Gasdermin-E protein (1-496AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A