

# DHFR Antibody, HRP conjugated



PACO26230

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Key enzyme in folate metabolism. Contributes to the de novo mitochondrial thymidylate biosynthesis pathway. Catalyzes an essential reaction for de novo glycine and purine synthesis, and for DNA precursor synthesis. Binds its own mRNA and that of DHFRL1.

**Gene ID:**

DHFR

**Uniprot**

P00374

**Synonyms:**

Dihydrofolate reductase (EC 1.5.1.3), DHFR

**Immunogen:**

Recombinant Human Dihydrofolate reductase protein (2-187AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A