

DNALI1 Antibody, Biotin conjugated



PACO26284

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a dynamic role in flagellar motility.

Gene ID:

DNALI1

Uniprot

O14645

Synonyms:

Axonemal dynein light intermediate polypeptide 1 (Inner dynein arm light chain, axonemal) (hp28), DNALI1

Immunogen:

Recombinant Human Axonemal dynein light intermediate polypeptide 1 protein (1-258AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A