

DNTTIP1 Antibody



PACO26305

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Shown to enhance TdT activity, in vitro.

Gene ID:

DNTTIP1

Uniprot

Q9H147

Synonyms:

Deoxynucleotidyltransferase terminal-interacting protein 1 (Terminal deoxynucleotidyltransferase-interacting factor 1) (Tdf1) (TdT-interacting factor 1), DNTTIP1, C20orf167 TDIF1


Immunogen:

Recombinant Human Deoxynucleotidyltransferase terminal-interacting protein 1 protein (1-329AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO26305 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO26305 at dilution of 1:100.