

EEF1D Antibody, HRP conjugated



PACO26382

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

EEF1D

Uniprot

P29692

Synonyms:

Elongation factor 1-delta (EF-1-delta) (Antigen NY-CO-4), EEF1D, EF1D

Immunogen:

Recombinant Human Elongation factor 1-delta protein (2-281AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A