

# EEF1D Antibody, FITC conjugated



PACO26383

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

EEF1D

**Uniprot**

P29692

**Synonyms:**

Elongation factor 1-delta (EF-1-delta) (Antigen NY-CO-4), EEF1D, EF1D

**Immunogen:**

Recombinant Human Elongation factor 1-delta protein (2-281AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A