

EEF1E1 Antibody, Biotin conjugated



PACO26388

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Positive modulator of ATM response to DNA damage.

Gene ID:

EEF1E1

Uniprot

O43324

Synonyms:

Eukaryotic translation elongation factor 1 epsilon-1 (Aminoacyl tRNA synthetase complex-interacting multifunctional protein 3) (Elongation factor p18) (Multisynthase complex auxiliary component p18), EEF1E1, AIMP3 P18

Immunogen:

Recombinant Human Eukaryotic translation elongation factor 1 epsilon-1 protein (2-174AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A