

# EMP3 Antibody, Biotin conjugated



PACO26472

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Probably involved in cell proliferation and cell-cell interactions.

**Gene ID:**

EMP3

**Uniprot**

P54852

**Synonyms:**

Epithelial membrane protein 3 (EMP-3) (Hematopoietic neural membrane protein 1) (HNMP-1) (Protein YMP), EMP3, YMP

**Immunogen:**

Recombinant Human Epithelial membrane protein 3 protein (25-65AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A