
Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:2000-1:10000,
IF:1:50-1:200

Protein Background:

Preferentially hydrolyzes nucleoside 5'-triphosphates. The order of activity with respect to possible substrates is UTP > GTP > CTP.

Gene ID:

ENTPD7

Uniprot

Q9NQZ7

Synonyms:

Ectonucleoside triphosphate diphosphohydrolase 7 (NTPDase 7) (EC 3.6.1. -)
(Lysosomal apyrase-like protein 1), ENTPD7, LALP1


Immunogen:

Recombinant Human Ectonucleoside triphosphate diphosphohydrolase 7 protein (51-350AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HeLa cells using PACO26493 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Western Blot

Positive WB detected in: HeLa whole cell lysate, Mouse heart tissue, Mouse spleen tissue, Mouse lung tissue, Mouse kidney tissue

All lanes: ENTPD7 antibody at 2.8µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 69 kDa

Observed band size: 69 kDa