

EPO Antibody, Biotin conjugated



PACO26519

Product Information

Size:

50ug

Reactivity:

Bovine

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Erythropoietin is the principal hormone involved in the regulation of erythrocyte differentiation and the maintenance of a physiological level of circulating erythrocyte mass.

Gene ID:

EPO

Uniprot

P48617

Synonyms:

Erythropoietin, EPO

Immunogen:

Recombinant Bovine Erythropoietin protein (26-192AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A