

ERG Antibody, FITC conjugated



PACO26531

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transcriptional regulator. May participate in transcriptional regulation through the recruitment of SETDB1 histone methyltransferase and subsequent modification of local chromatin structure.

Gene ID:

ERG

Uniprot

P11308

Synonyms:

Transcriptional regulator ERG (Transforming protein ERG), ERG

Immunogen:

Recombinant Human Transcriptional regulator ERG protein (1-486AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A