

# ERG Antibody, Biotin conjugated



PACO26532

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Transcriptional regulator. May participate in transcriptional regulation through the recruitment of SETDB1 histone methyltransferase and subsequent modification of local chromatin structure.

**Gene ID:**

ERG

**Uniprot**

P11308

**Synonyms:**

Transcriptional regulator ERG (Transforming protein ERG), ERG

**Immunogen:**

Recombinant Human Transcriptional regulator ERG protein (1-486AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A