

ETFDH Antibody



PACO26561

Product Information

Size:

50ug

Reactivity:

Pig

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Accepts electrons from ETF and reduces ubiquinone.

Gene ID:

ETFDH

Uniprot

P55931

Synonyms:

Electron transfer flavoprotein-ubiquinone oxidoreductase, mitochondrial (ETF-QO) (ETF-ubiquinone oxidoreductase) (EC 1.5.5.1) (Electron-transferring-flavoprotein dehydrogenase) (ETF dehydrogenase), ETFDH

Immunogen:

Recombinant Pig Electron transfer flavoprotein-ubiquinone oxidoreductase, mitochondrial protein (34-617AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A