

EVI2A Antibody



PACO26577

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

May complex with itself or/and other proteins within the membrane, to function as part of a cell-surface receptor.

Gene ID:

EVI2A

Uniprot

P22794

Synonyms:

Protein EVI2A (Ecotropic viral integration site 2A protein homolog) (EVI-2A), EVI2A, EVDA EVI2


Immunogen:

Recombinant Human Protein EVI2A protein (31-133AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO26577 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO26577 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).