

# F8 Antibody, HRP conjugated



PACO26618

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Factor VIII, along with calcium and phospholipid, acts as a cofactor for factor IXa when it converts factor X to the activated form, factor Xa.

**Gene ID:**

F8

**Uniprot**

P00451

**Synonyms:**

Coagulation factor VIII (Antihemophilic factor) (AHF) (Procoagulant component)  
[Cleaved into: Factor VIIIa heavy chain, 200 kDa isoform; Factor VIIIa heavy chain, 92 kDa isoform; Factor VIII B chain; Factor VIIIa light chain], F8, F8C

**Immunogen:**

Recombinant Human Coagulation factor VIII protein (1-216AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A