

TAGLN Antibody, Biotin conjugated



PACO26656

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Actin cross-linking/gelling protein By similarity. Involved in calcium interactions and contractile properties of the cell that may contribute to replicative senescence.

Gene ID:

TAGLN

Uniprot

Q01995

Synonyms:

Transgelin (22 kDa actin-binding protein) (Protein WS3-10) (Smooth muscle protein 22-alpha) (SM22-alpha), TAGLN, SM22 WS3-10

Immunogen:

Recombinant Human Transgelin protein (2-201AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A