

FAM3B Antibody, FITC conjugated



PACO26675

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Induces apoptosis of alpha and beta cells in a dose- and time-dependent manner.

Gene ID:

FAM3B

Uniprot

P58499

Synonyms:

Protein FAM3B (Cytokine-like protein 2-21) (Pancreatic-derived factor) (PANDER), FAM3B, C21orf11 C21orf76

Immunogen:

Recombinant Human Protein FAM3B protein (31-235AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A