

FAM92A Antibody, HRP conjugated



PACO26678

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	FAM92A
Source:	Uniprot
Rabbit	A1XBS5
Isotype:	Synonyms:
IgG	Protein FAM92A, FAM92A, FAM92A1
Applications:	Immunogen:
ELISA	Recombinant Human Protein FAM92A protein (1-271AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A