

CDA Antibody, Biotin conjugated



PACO26696

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

This enzyme scavenges exogenous and endogenous cytidine.

Gene ID:

CDA

Uniprot

P32320

Synonyms:

Cytidine deaminase (EC 3.5.4.5) (Cytidine aminohydrolase), CDA, CDD

Immunogen:

Recombinant Human Cytidine deaminase protein (1-146AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A