

FDFT1 Antibody, Biotin conjugated



PACO26720

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

endoplasmic reticulum, endoplasmic reticulum membrane, integral component of membrane, farnesyl-diphosphate farnesyltransferase activity, cellular lipid metabolic process, cholesterol biosynthetic process, small molecule metabolic process, steroid biosynthetic process.

Gene ID:

FDFT1

Uniprot

P37268

Synonyms:

Squalene synthase (SQS) (SS) (EC 2.5.1.21) (FPP: FPP farnesyltransferase) (Farnesyl-diphosphate farnesyltransferase), FDFT1

Immunogen:

Recombinant Human Squalene synthase protein (2-417AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A