

FDXR Antibody, FITC conjugated



PACO26723

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Serves as the first electron transfer protein in all the mitochondrial P450 systems. Including cholesterol side chain cleavage in all steroidogenic tissues, steroid 11-beta hydroxylation in the adrenal cortex, 25-OH-vitamin D3-24 hydroxylation in the kidney, and sterol C-27 hydroxylation in the liver.

Gene ID:

FDXR

Uniprot

P22570

Synonyms:

NADPH: adrenodoxin oxidoreductase, mitochondrial (AR) (Adrenodoxin reductase) (EC 1.18.1.6) (Ferredoxin--NADP(+) reductase) (Ferredoxin reductase), FDXR, ADXR

Immunogen:

Recombinant Human NADPH: adrenodoxin oxidoreductase, mitochondrial protein (33-451AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A