

FFAR2 Antibody, Biotin conjugated



PACO26728

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

FFAR2

Uniprot

O15552

Synonyms:

Free fatty acid, receptor 2 (G-protein coupled receptor 43), FFAR2, FFA2 GPCR43 GPR43

Immunogen:

Recombinant Human Free fatty acid, receptor 2 protein (277-330AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A