

Fgf1 Antibody, FITC conjugated



PACO26743

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays an important role in the regulation of cell survival, cell division, angiogenesis, cell differentiation and cell migration. Functions as potent mitogen in vitro.

Gene ID:

Fgf1

Uniprot

P61148

Synonyms:

Fibroblast growth factor 1 (FGF-1) (acid, c fibroblast growth factor) (aFGF) (Heparin-binding growth factor 1) (HBGF-1), Fgf1, Fgf-1 Fgfa

Immunogen:

Recombinant Mouse Fibroblast growth factor 1 protein (16-155AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A