

Fgf21 Antibody, Biotin conjugated



PACO26752

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Stimulates glucose uptake in differentiated adipocytes via the induction of glucose transporter SLC2A1/GLUT1 expression (but not SLC2A4/GLUT4 expression). Activity probably requires the presence of KLB.

Gene ID:

Fgf21

Uniprot

Q9JJN1

Synonyms:

Fibroblast growth factor 21 (FGF-21), Fgf21

Immunogen:

Recombinant Mouse Fibroblast growth factor 21 protein (29-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A