

FGFR1OP2 Antibody



PACO26761

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

May be involved in wound healing pathway.

Gene ID:

FGFR1OP2

Uniprot

Q9NVK5

Synonyms:

FGFR1 oncogene partner 2, FGFR1OP2

Immunogen:

Recombinant Human FGFR1 oncogene partner 2 protein (1-172AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: FGFR1OP2 antibody at 6 μ g/ml

Lane 1: Mouse gonadal tissue

Lane 2: Human placenta tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 30, 25, 21 kDa

Observed band size: 30 kDa

.