

FGFR1OP2 Antibody, FITC conjugated



PACO26763

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in wound healing pathway.

Gene ID:

FGFR1OP2

Uniprot

Q9NVK5

Synonyms:

FGFR1 oncogene partner 2, FGFR1OP2

Immunogen:

Recombinant Human FGFR1 oncogene partner 2 protein (1-172AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A