

# LTF Antibody, HRP conjugated



PACO26786

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Transferrins are iron binding transport proteins which can bind two Fe<sup>3+</sup> ions in association with the binding of an anion, usually bicarbonate. Ref.23

Ref.30 Lactotransferrin has antimicrobial activity which depends on the extracellular cation concentration.

**Gene ID:**

LTF

**Uniprot**

P02788

**Synonyms:**

Lactotransferrin (Lactoferrin) (EC 3.4.21. -) (Growth-inhibiting protein 12) (Talalactoferrin) [Cleaved into: Lactoferricin-H (Lfcin-H); Kaliocin-1; Lactoferroxin-A; Lactoferroxin-B; Lactoferroxin-C], LTF, GIG12 LF

**Immunogen:**

Human Lactotransferrin protein (Native Protein).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A