

LTF Antibody, FITC conjugated



PACO26787

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transferrins are iron binding transport proteins which can bind two Fe³⁺ ions in association with the binding of an anion, usually bicarbonate. Ref.23

Ref.30 Lactotransferrin has antimicrobial activity which depends on the extracellular cation concentration.

Gene ID:

LTF

Uniprot

P02788

Synonyms:

Lactotransferrin (Lactoferrin) (EC 3.4.21. -) (Growth-inhibiting protein 12) (Talalactoferrin) [Cleaved into: Lactoferricin-H (Lfcin-H); Kaliocin-1; Lactoferroxin-A; Lactoferroxin-B; Lactoferroxin-C], LTF, GIG12 LF

Immunogen:

Human Lactotransferrin protein (Native Protein).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A