

FLG Antibody, Biotin conjugated



PACO26792

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Aggregates keratin intermediate filaments and promotes disulfide-bond formation among the intermediate filaments during terminal differentiation of mammalian epidermis.

Gene ID:

FLG

Uniprot

P20930

Synonyms:

Filaggrin, FLG

Immunogen:

Recombinant Human Filaggrin protein (3838-4061AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A